

# Handling Belt

Code : MF / HB1



This soft padded patient handling belt can be adjusted to fit around the waist of your patient and provide you with a choice of secure hand holds while performing assisted manoeuvres or transfers.

- Designed to enable better posture for staff while giving more control over the patient.
- Ergonomically designed the belt is soft and padded making it comfortable to wear.
- Multiple handles on the belt provide the best possible grip in a convenient position which helps caregivers maintain a safe posture.
- non-slip inner lining to help prevent the belt sliding upwards during a transfer (sit to stand, or assisted walking).
- Adjustment in diameter is achieved by simply extending the webbing girdle to suit small to large patients.
- Strong Handles are comfortably shaped, and placed at oblique angles to aid comfort and enhance control for the carer.
- Made from washable materials making it easy to keep hygienically clean.
- Can be cleaned with hot soapy water and rinse.

## Used in typical situations such as :

Sit to stand, help when walk training, transferring on and off the toilet, in and out of the bath, lounge, wheel chair or vehicle. In fact anywhere a care giver is performing assisted manoeuvres or transfers this handling belt will make your task easier and safer.

Handling belts are popular with hospitals, stroke units, general ward use, ambulance staff, rescue, community based staff, care at home, Institutional Care, nursing and residential homes.

## Note :

Care should be taken not to use this belt for lifting and should never be worn by the care giver for the patient to hold onto as this could result in severe injury to the care giver should the patient fall.



## Specifications

Size of handles	90mm wide x 25mm deep
Number of handles	4
Belt diameter adjustment	700 mm - 1340 mm
Width of padded belt	190 mm - 130 mm
Width of adjustable strap	50 mm
Max user	120kg with larger sizes on

## Multifit

Hospital Supplies Ltd

Multifit Hospital Supplies Ltd

Ph 0800 102 803

Fax 0800 102 800

info@multifit.co.nz

www.multifit.co.nz